

## Spectrometer (S.S Scale)



**Product Categories:** [Microscope](#), [Microscopes](#)

**Product Tags:** [biology lab equipment](#)

**Product Page:** <https://www.labappara.com/product/spectrometer-s-s-scale/>

### Product Description

Spectrometer (S.S Scale) - Laboratory Equipment

Deluxe with Graticule, specification as above in wooden box.

For Introductory class work

- A Basic measuring instrument for quantitative spectroscopy experiments, this instrument is mounted on a stable cast iron base with attached collimator and a rotating telescope platform and graduated circle.
- A 2.4 cm. Diameter achromatic optical system with a 10 x Ramsden eyepiece & cross line graticule.
- An 8 cm. Diameter prism/grating Table, fully adjustable with prism & grating holders and engraved ring pattern. Protected Circle 15.0 cm with

**Stainless steel Scale**