

## Spirometer



**Product Categories:** [Pathology Lab](#)

**Product Tags:** [Laboratory Equipment](#)

**Product Page:** <https://www.labappara.com/product/spirometer/>

### Product Description

#### **Spirometer- Laboratory Equipment**

6 litre capacity for vital capacity, or use as a gasometer with chain compensated counter balance to the float balanced through cut its movement.

Pulley calibrated to denote volume inlet and outlet tubes, side cock for gases, cock for draining off water strong construction. Complete with corrugated tube mouth piece & recording lever.

An instrument for measuring the air capacity of the lungs.

#### **Feature -**

6 Liter capacities with chain compensate counter balance to float thought pulley.

Calibrate to donate volume, inlet and outlet tubes.

Spirometer is made consists of stainless steel chamber.