

## Storage Vial-Internal Thread



**Product Categories:** [Laboratory Plastic Ware](#), [Vials](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#)

**Product Page:** <https://www.labappara.com/product/storage-vial-internal-thread/>

### Product Description

Storage Vial-Internal Thread - Scientific Laboratory Equipment

Another addition to the range of storage vials is this vial which has internal thread in it. The cap fitted with an O-ring screws in to the vial and makes it a leakproof storage vial.

Material : Polypropylene

### Particulars

Packing

2.0 ml

500 Pcs