


Study Of OPTO - Coupler



Product Categories: [Characteristics And Application](#), [Characteristics And Application](#), [Educational equipments](#), [Electronics](#), [Engineering Equipment](#)

Product Tags: [electronics](#)

Product Page: <https://www.labappara.com/product/study-opto-coupler/>

Product Description

Study Of OPTO - Coupler - Electronics Laboratory Instrument

Study Of OPTO - Coupler with build in two regulated power supply.

And with four meters for measurement of voltage & current.

Complete in all respect.

Objective:

To construct opto coupler circuit

To study its voltage & its current response

To study its frequency response

To plot its characteristics curve

Technical Specification : Built in power supply : +5V & +15V DC

Mains : 230V AC

Dimension : 27 CMS X 17 CMS [metal cabinet]

Weight : 2KGS APPX