

Test Tube Holder (Iron)



Product Categories: [Chemistry Equipments](#), [General Chemistry Labware](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/test-tube-holder-iron/>

Product Description

Test Tube Holder (Iron) - Chemistry Laboratory Equipment

A test tube holder is used to hold test tubes. It is used for holding a test tube in place when the tube is hot or should not be touched. It's maximum opening is 20 to 25 mm. It's length is 200mm, full accuracy. It can also be used to hold the test tube when the test tube is being heated.

Features:

A test tube holder is an apparatus designed to support and store test tubes during a scientific experiment.

Test tube holders may also be referred to as test tube racks.

Thick flat steel wire, chrome plated have thumb sliding device for opening and closing, with polished wooden handle.

Supreme quality

Better grip for tubes

Durability

AVAILABILITY (ADD DIFF PICS)

Flat form with wooden handle.

Steel zinc plated.

Cross pattern iron Z.P