

TEST TUBES



Product Categories: [Laboratory Glassware](#), [Test Tubes / Culture Tubes](#)

Product Tags: [chemical glassware](#), [chemistry lab instruments](#), [scientific glassware](#)

Product Page: <https://www.labappara.com/product/test-tubes/>

Product Description

Test Tubes - Laboratory Glassware Equipment

A test tube, also known as a culture tube or sample tube, is a common piece of laboratory glass ware consisting of a finger-like length of glass Test tubes are available in a multitude of lengths and widths Test tubes are available in a multitude of lengths and widths The top often features a flared lip to aid pouring out the contents;

Available In

Borosilicate glass tubing

Round bottom

With rim / Rimless

Plain.

Graduation in Amber color.

Made from Neutral Hard Glass.

Capacity O.D. x Length Minimum

ml (appx) Size Qty/Pack

- 3 10x75 100
- 5 12x75 100
- 7 12x100 100
- 10 15x125 100
- 15 15x150 100
- 25 18x150 100
- 35 25x100 50
- 50 25x150 50
- 75 25x200 50
- 100 32x200 25
- 175 38x200 25