

Thermistor Characteristics Apparatus



Product Categories: [Characteristics And Application](#), [Characteristics And Application](#), [Electronics](#), [Engineering Equipment](#)

Product Page: <https://www.labappara.com/product/thermistor-characteristics-apparatus/>

Product Description

Thermistor Characteristics Apparatus

Features :

Instrument comprises of DC Regulated Power Supply 0-10 VDC/150mA, two round meters for voltage & current measurement, one thermistor is kept in Oven. Connections of Supplies, Meters & Thermistor brought out at 4 mm Sockets, Also provided with Oven with temperature control knob.