

## Transistor Characteristic Apparatus



**Product Categories:** [Characteristics And Application](#), [Educational equipments](#), [Electronics](#)

**Product Tags:** [electronics](#)

**Product Page:**

<https://www.labappara.com/product/transistor-characteristic-apparatus/>

### Product Description

Transistor Characteristic Apparatus - Electronics Laboratory Instruments  
Transistor Characteristic Apparatus with build in stabilized power supply.  
Workable on 230 volts.

AC fitted with PNP transistor or AC 128 on front panel.

Different electrodes of transistor & power supply.

Output and meter connections are brought in panel with separate terminals.

With in four meters

### **Specification :**

Instrument comprises of Two DC Regulated Power Supplies 0-10 V DC/150mA 0-1 V DC/150 mA, four acrylic meters for voltage & current measurement , one PNP & one NPN Transistor mounted on the panel, connections of Supplies , Transistors brought out at 4 mm Sockets. NPN or NPN transistor can be selected by use of toggle switch. Special design of double tapered powder coated box make it

very Elegant, attractive & it also gives proper viewing of analog meters. Set of 10 Single point interconnection leads & comprehensive User Manual is provided.

**Features:**

Rugged construction

Compact design

Easy to use

Highly accurate