

Triode Valve Characteristics Apparatus



Product Categories: [E/M Apparatus & Valve Experiments](#), [E/M Apparatus & Valve Experiments](#), [Electronics](#), [Engineering Equipment](#)

Product Page:

<https://www.labappara.com/product/triode-valve-characteristics-apparatus/>

Product Description

Triode Valve Characteristics Apparatus

Features :

Instrument comprises of DC Regulated Power Supply 0-300VDC/ 30mA & 0-25VDC/10mA, AC Power Supply 6.3 VAC/1Amp, two round meters for voltages & current measurement, valve diagram is printed on front panel, connections of supplies & valve brought out at 4mm socket.