

## Triple Beam Balance



**Product Categories:** [Balance & Scale](#), [Chemistry Equipments](#)

**Product Tags:** [chemistry lab instruments](#), [physics lab instruments](#), [science lab equipment](#)

**Product Page:** <https://www.labappara.com/product/triple-beam-balance/>

### Product Description

Triple Beam Balance - Laboratory Instrument

The Triple beam balance is a rugged, inexpensive weighing instrument.

It is named for the three beams which carry weights. The middle beam reads only in 100 g increments.

Triple beam balance is an instrument used to measure mass very precisely.

The name refers to the three beams including the middle beam which is the largest size,

The front beam which is the medium size,

And the far beam which is the smallest size.

**Capacity - 111 to 310 gm Feature -**

For direct reading

With sliding riders

With sunmica base.

Improved model

Export Quality