

# TTL-CMOS, CMOS-TTL Interfacing



**Product Categories:** [Digital Electronics Lab](#), [Engineering Equipment](#)

**Product Page:** <https://www.labappara.com/product/ttl-cmos-cmos-ttl-interfacing/>

## Product Description

TTL-CMOS, CMOS-TTL Interfacing

**Features :**

Instrument comprises of Two DC Regulated Power Supplies of 5V/150mA  
12V/150mA Output indicator for CMOS & TTL individually.

Built in pull up network of resistor.

ICs are mounted on front panel & appropriate test point are provided on front panel.

One CMOS & One TTL IC are provided for experiment.