

## V.D.R.L Shaker



**Product Categories:** [Pathology Lab](#)

**Product Tags:** [Laboratory Equipment](#)

**Product Page:** <https://www.labappara.com/product/v-d-r-l-shaker/>

### Product Description

#### V.D.R.L Shaker - Laboratory Equipment

VDRL shaker is an essential equipment for a instruments laboratory, where it is used for the mixing of syphilis testing card, emulsion reaction, and complement fixation reaction while can also be used to blends of clinical diagnostic tests requiring gentle agitation.

Shakers are Very useful for V.D.R.L typing grouping as well another agglutination tests.

V.D.R.L Shaker is also for mixing of solution in small flasks, beakers, bottles or similar containers.

Shakers are operated on 230 vol AC.

Compact shaker for mixing blood samples.

Timer options for accurate mixing time.

Digital RPM controller and digital timer available for exact results.

Stable at high speeds.

Wide speed range.

Separate platform attachment available for mixing more than one samples.

**Technical Specification**

Type of movement	Rotational
Speed range	0- 350 RPM
Operational mode	Continuous operation
Timer option	Yes
Dimensions	37 x 37 x 18
Weight	11 kgs
Permissible ambient temperature	5 - 40 deg.
Permissible relative moisture	80%
Motor type	PMDC motor
Maximum load	6- kg
Speed display	None
Supply	220 V, 50 Hz.