

Vibrating Screen



Product Categories: [Engineering Equipment](#), [Mechanical Operation Lab](#)

Product Page: <https://www.labappara.com/product/vibrating-screen/>

Product Description

Vibrating Screen

Vibrating Screen are those, which rapidly vibrated with small amplitude, are less likely to blind than gyrating screens. The vibrations may be generated mechanically or electrically. Mechanical vibrations are usually transmitted from high-speed eccentrics to the casing of the unit and from there to inclined screens. The setup consists of three different screen with 380 mm width and 600 mm length fitted on a rigid M.S. Frame, feed hopper, eccentric etc. Three screens of different sizes are fitted one above the other. The drives are given from an electric motor.

Demonstration of vibrating screen.

A good quality painted rigid ms structure is provided to support all the parts

FEATURES :

Simple and dependable Eccentric type system.

High screening capacity.

High durability.

No transmission of screen panels.

Ready changing of screen plates.

Rigid and vibrating resistance screen body.

