

Voltage Stabilization Characteristics Of Zener Diode Apparatus



Product Categories: [Characteristics And Application](#), [Electronics](#)

Product Page:

<https://www.labappara.com/product/voltage-stabilization-characteristics-of-zener-diode-apparatus/>

Product Description

Voltage Stabilization Characteristics Of Zener Diode Apparatus

Features :

Instrument comprises of One DC regulated Power Supply, two round meters to measure input voltage & output voltage, one Zener Diode mounted behind the panel, built in load resistance selectable using band switch & Connections brought out at terminals/sockets.