

## Voltage Stabilization Characteristics Of Zener Diode Apparatus



**Product Categories:** [Characteristics And Application](#), [Electronics](#)

**Product Page:**

<https://www.labappara.com/product/voltage-stabilization-characteristics-of-zener-diode-apparatus/>

### Product Description

#### Voltage Stabilization Characteristics Of Zener Diode Apparatus

**Features :**

Instrument comprises of One DC regulated Power Supply, two round meters to measure input voltage & output voltage, one Zener Diode mounted behind the panel, built in load resistance selectable using band switch & Connections brought out at terminals/sockets.