

Voltmeter



Product Categories: [Lab Meters](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/voltmeter/>

Product Description

VOLTMETER- Physics Lab Equipments

A voltmeter, also known as a voltage meter, is an instrument used for measuring the potential difference, or voltage, between two points in an electrical or electronic circuit. It is generally used in all types of Laboratories such as School, college, university, Research etc.

Dc mounted in beautiful Bakelite / acrylic stand with accuracy $\pm 0\%$ F.S.D.

With lock Type terminal packed in card box.

All common ranges available.

Available Models – MO65, MR80, MR100