

## Wash Bottle (polyethylene)



**Product Categories:** [Bottles And Assesories](#), [Laboratory Plastic Ware](#)

**Product Tags:** [Laboratory Equipments](#), [Laboratory plastic ware](#)

**Product Page:** <https://www.labappara.com/product/wash-bottle-polyethylene/>

### Product Description

Wash Bottle (polyethylene) - Scientific Laboratory Equipment

Poly-lab wash bottles are made of low density polyethylene which gives these bottles a translucent & unbreakable character. Solutions can be poured easily by squeezing the bottle lightly. The cap is fitted with a flexible polyethylene delivery tube which can be aimed wherever required. After squeezing the Wash Bottle lightly solutions can be easily poured.

Material : Polyethylene

Particulars	Packing
125 ml	12 Pcs
250 ml	12 Pcs
500 ml	6 Pcs
1000 ml	6 Pcs