

## Wash bottle



**Product Categories:** [Bottles And Asssesories](#), [Laboratory Plastic Ware](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#)

**Product Page:** <https://www.labappara.com/product/wash-bottle/>

### Product Description

Wash bottle - Scientific Laboratory Equipment

Poly-lab wash bottles are made of low density polyethylene which gives these bottles a translucent & unbreakable character. Solutions can be poured easily by squeezing the bottle lightly. The cap is fitted with a flexible polyethylene delivery tube which can be aimed wherever required.

Wash Bottle are Translucent and unbreakable character. After squeezing the Wash Bottle lightly solutions can be easily poured. The Cap is fitted with a Flexible Polyethylene Delivery Tube which can be aimed wherever required.

Material : Polyethylene