

Watch Glass



Product Categories: [Chemistry Equipments](#), [General Chemistry Labware](#), [Optical Instruments](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#), [science lab equipment](#)

Product Page: <https://www.labappara.com/product/watch-glass/>

Product Description

Watch Glass - Scientific Equipment

Watch glass is used for multiple purposes in laboratories. Watch glass is a circular concave piece of glass used in chemistry as a surface to evaporate a liquid, to hold solids while being weighed, for heating a small amount of substance and as a cover for a beaker.

Available in size – 50 mm, 60 mm, 75 mm, 100 mm, 125 mm and 150 mm. Diameters.

Feature -

Uniformity module ground edge.

Well packed in card box.

Watch Glasses Polished, ground edges, made from selected glass sheet.

Crack resistance

Perfect finish

Highly durable