

Wide Field Eyepiece



Product Categories: [Eye Piece](#), [Microscope](#), [Microscope Assecories](#), [Optical Instruments](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/wide-field-eye-piece/>

Product Description

Wide Field Eyepiece - Scientific Equipment

Anti-reflection, anti-fungus coated 10 x wide-field pair compatible with most of all microscopes.

With cross line reticulated suitable for spectrometer.

Consists of two Plano-convex lenses positioned with both convex faces toward the object.

Designed for use with low power achromatic objectives.

Our entire line of eyepieces will fit any international standard microscope with a 23 mm inner diameter eyepiece tube.