

## Windmill kit (wooden)



**Product Categories:** [Educational equipments](#), [Geography & Working Models](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#)

**Product Page:** <https://www.labappara.com/product/wind-mill-kit-wooden/>

### Product Description

#### Windmill kit (wooden) - Scientific Laboratory Educational Equipment

This is made wooden kit to demonstrate the generation of electric energy by a wind energy (kinetic energy) or generation of electricity by moving a coil in a magnet.

When air strikes the generator blade, it turns the coil of small generator to produces electricity and lights up a LED and a connected voltmeter indicates the amount of current being produced by the generator