


## Wire Gauge



**Product Categories:** [Chemistry Equipments](#), [Educational equipments](#), [General Chemistry Labware](#), [Physics Equipments](#)

**Product Tags:** [Chemistry Equipments](#)

**Product Page:** <https://www.labappara.com/product/wire-gauge/>

### Product Description

Wire Gauge - Chemistry Equipment

**Wire gauge** is most commonly used in chemistry Laboratories. It is put on top of the Bunsen Burner.

Then then some glassware like a breaker is put on the Bunsen burner to heat.

Material: Made from G.I. wire with Wire Gauge With asbestos center, extra strong quality in squares.

**Features:**

Heat resistance

Easy to fit

Safe usage