


## Zener Diode Characteristic Apparatus



**Product Categories:** [Characteristics And Application](#), [Educational equipments](#), [Electronics](#)

**Product Tags:** [electronics](#)

**Product Page:**

<https://www.labappara.com/product/zener-diode-characteristic-apparatus/>

### Product Description

Zener Diode Characteristic Apparatus - Electronics Equipment

Zener Diode Characteristic Apparatus with build in stabilized power supply.

Workable on 230 volts , AC fitted with zener diode 8.2 volts on front panel connections of zener diode.

Power supply , output and meters connections are brought out on the panel.

Output 0-10 volts DC .

Continuously variable.

### Specification :

Instrument comprises of DC Regulated Power supply 0-15 V DC 150 mA,

Two round meters for voltage & current measurement,

3 Zener Diodes mounted behind the panel, connections of Supplies,

Meters & zener Diode brought out at 4 mm Sockets.

### Features:

Easy maintenance

Accuracy

Robustness

Sturdy design

Long working life

Perfect functioning